MARYLAND STATE DEPARTMENT OF HEALTH AND MENTAL HYGIENE LABORATORIES ADMINISTRATION

PLEASE READ - REVISED INSTRUCTIONS

USING THE SEROLOGICAL TESTING FORM TO REQUEST HIV TEST

NOTE: References to the CTR site# and CTR # ID sticker apply ONLY to participants in the statewide CTR System

I. DIRECTIONS FOR COLLECING BLOOD AND LABELING THE SPECIMEN TUBE

Gloves should always be worn when collecting blood from individuals. Use a disposable Vacutainer system to collect blood. Collect 5-7 ml of blood in a red stopper Vacutainer tube and allow blood to clot. All specimen tubes should be labeled with the patient name. Please affix patient ID labels to the specimen tube.

IF THE VACUTAINER ADAPTOR AND NEEDLE ARE NOT AVAILABLE FOR COLLECTION OF BLOOD: Use only disposable syringes and needles to collect blood. Transfer blood to Vacutainer tube by inserting the needle through the stopper and allowing the tube vacuum to pull the blood from the syringe into the tube.

DO NOT REMOVE THE STOPPER FROM THE VACUTAINER TUBE OR FORCE BLOOD INTO THE TUBE. If the stopper is removed, it cannot be replaced securely for transport. DO NOT ATTEMPT TO REMOVE THE NEEDLE FROM THE SYRINGE OR BEND/CUT THE NEEDLE PRIOR TO DISPOSAL.

II. FILLING OUT LAB REQUEST SLIP

Please provide all available information on the lab request form. For submitter information refer to Comment Boxes A and B on the sample form. For patient information refer to Comment Box C. If a CTR ID number sticker was used on the specimen tube then place the same label on the lab request form, refer to Comment Box D. Include both the date and time of specimen collection on lab request form; refer to Comment Box E. Indicate if this is an initial test or follow-up test by marking 1st specimen or 2nd specimen respectively. Include previous test information in designated area of the form. CTR sites should provide the site #; refer to Comment Box F. The last copy of the Serological Testing Form should be retained at the collection site.

Use printed patient ID labels whenever possible on both the lab request form and specimen tube to ensure the spelling of the patient name and patient ID numbers match exactly.

HIV SECTION OF THE SEROLOGICAL TESTING FORM (See red outlined area and refer to Comment Box G)

III. SUBMITTING SPECIMENS

The DHMH Laboratories Administration has approved the use of bio-bags for the transport of HIV specimens via courier when the blood has been collected in a <u>plastic</u> vacutainer tube. **DO NOT** transport **glass** vacutainer tubes in bio-bags. Place specimen tube into the inner pocket of the bio-bag and securely seal the bag, then place the lab request form in the outer pocket of the bio-bag.

ONLY ONE SPECIMEN TUBE SHOULD BE SUBMITTED PER SINGLE BIO-BAG. DO NOT BUNDLE, SECURE WITH RUBBER BANDS, TAPE OR STAPLE THE BIO-BAGS.

When transporting through US Mail, UPS or FEDEX, federally approved Biohazard packaging must be used (refer to U.S. Postal Service Regulations, Department of Transportation Regulations and International Air Transport Association Regulations).

Acceptance criteria for HIV Testing: Specimens must be received within 3 days if stored and shipped at room temperature or 7 days if stored refrigerated and shipped on cold packs/ in cooler. Specimens should be stored refrigerated whenever possible and submitted for testing immediately after collection to ensure a satisfactory test specimen. Record specimen storage and transport conditions on the Serological Testing Form in the designated HIV Section outlined in red and refer to Comment Box G. A separate specimen must be submitted for HIV testing. SPECIMENS GREATER THAN 7 DAYS ARE UNACCEPTABLE FOR TESTING.

STATE LAB Use Only

A. Complete fields of the 1205 443-submitter information section or place preprinted submitter information label here.

1205

boratories/

MARYLAND

C. Complete fields of the patient information section or place preprinted patient ID label here.

	□EH □FP□MTY/PN□NOD □STD □TB □CD □COR			
TYPE OR PRINT REQUIRED INFORMATION OR PLACE LABELS ON ALL THREE COPIES	Health Care Provider		Patient SS# (last 4 digits):	
	Address		Last Name	□ SR □ JR □ Other
RMA	City County		First Name	M.I.
FOF	State Zip Code		Date of Birth (mm/dd/yyyy)	1 1
≥ E	Contact Name:		Address	
ALI	Phone# Fax#		City	County
O S	rest reducet rationized by:		TRAB box or include on	Zip Code
TRE	Sex: Male Female Transgender M to F your label/sta		amp the name/credential	
AB			rized provider ordering the	er Pacific Islander 🖂 White
RP	MRN/Case # DOC # test.			_ab#
LA O	Date Collected:			story:
TY P	Previous Test Done? no yes Name of Test Date 1st 2nd 3rd State Lab Number:			
	Name of Test Date			
	Onset Date: Clinical Illness/Sympt			toms:
♣ SPE	CIMEN SOURCE CODE	♣ SPECIMEN S		♣ SPECIMEN SOURCE CODE
	Arbovir E. Collection date and	Herpes Si	mplex Virus (HSV) Types 1&2	♣ LAVENDER TOP TUBE REQUIRED
Mandatory: Onset time must be entered.		Legionella	l	Hemoglobin Disorders
Arbovirus Endemic Panel (WNV, EEE, SLE, LAC)		Leptospira	3	Blood transfusion? (last 4 months)
	Arbovirus Travel-Associated Panel	Lyme Disease		□ yes □ no
(Chikungunya, Dengue)		*MMRV Immunity Screen: [Measles		Prenatal screen? ☐ yes ☐ no
	sed on information provided PCR and/or nunological assays will be performed.	(Rubeola), Mumps, Rubella, Varicella		Father of baby screen? ☐ yes ☐ no
555	quired information, check all that apply:	(Chickenpox) IgG Ab only]		Guardian's name if patient is a minor:
DIAGNOSIS: □Aseptic Meningitis □Encephalitis □other □ SYMPTOMS: □headache □fever □stiff neck □altered mental state □muscle weakness □rash □other □ ILLNESS FATAL? □ yes □no TRAVEL HISTORY (dates and places)		Mononucleosis - Infectious		Name of mother of "at risk" baby:
		*Mumps Immunity Screen		1
		Mycoplasma		
		Rocky Mountain Spotted Fever (RMSF)		RESTRICTED TEST Pre-Approved Submitters Only Submit a separate specimen for HIV Instructions go to:
		*Rabies (RFFIT) (*List vaccination dates above)		
		*Rubella Immunity Screen		
		*Rubeola ((Measles) Immunity Screen	http://dhmh.maryland.gov/laboratories/
IMMUNIZATIONS: Yellow fever? ☐ yes ☐ no Flavivirus? ☐ yes ☐ no		Schistosor	na	B HIN
		Strongyloides		Country of Origin
IMI	MUNOCOMPROMISED? ☐ yes ☐no	Syphilis - F	Previously treated? ☐ yes ☐ no	Rapid Test: ☐ Reactive ☐ Negative
		Toxoplasm	na /	Date:
	Aspergillus	Tularemia		Specimen stored refrigerated (2°-8°c)
	chlamydia (group antigen IgG) G. Must submit separate specimen and write / after collection.			
-	Cryptococcal (antigen) specimen code (B or S) next to HIV test Specimen transported on cold packs			
+	Cytomegalovirus (CMV) order. **Note X or check marks are not			□ yes □ no
	Ehrlichia acceptable.**			SPECIMEN SOURCE CODE:
	Epstein-Barr Virus (EBV)		angements have been made with the	PLACE CODE IN BOX NEXT TO TEST
Hepatitis A Screen (IgM Ab only, acute infection)		following DHMH Labs Administration employee:		B Blood (5 ml) CSF Cerebrospinal Fluid
	Call lab (443-681-3889) prior to submitting	ciripioyee	2.	L Lavender Top Tube
	Hepatitis B Screen (HBs antigen only)			P Plasma
	Prenatal patient? ☐ yes ☐ no			S Serum (1 ml per test)
	Hepatitis B Panel: (HBsAg, HBsAb)	F. Writ	e CTR Site# here.	UR Urine
	Hepatitis B post vaccine (HBsAb)	2		D. Place CTR ID# Sticker
	Hepatitis C screen (HCV Ab only)		accination History above*	here.
DHMH 4677 Revised 08/15 Original				

CLINIC CODES

EH - Employee Health

FP - Family Planning

MTY/PN - Maternity/Prenatal

NOD - Nurse of Day

STD/STI - Sexually Transmitted Disease/Infections

CD - Communicable Disease

COR - Correctional Facility

Do not mark a box if clinic type does not apply

COMPLETING FORM

Type or print legibly

Printed labels are recommended

Place printed labels on all copies of form

Write the person's name that is authorized to order test in the box provided

Press firmly - three part form

Collection date & time are required by Law

WRITE SPECIMEN CODE in box next to test

*Specimens/samples can not be processed without a requested test.

VACCINATION HISTORY

List vaccination dates for all Rabies, Hepatitis B and MMRV (Mumps, Measles, Rubella, and Varicella) test request.

Rabies Vaccination history is required for all RFFIT test requests.

HIV TESTING

Include previous HIV Test information in the top section under Previous Test done Submit a separate specimen for HIV testing when multiple tests are ordered on the one form

Questions/comments on the use of the specimen bags/storage/shipping or completing the form contact:

Accessioning Unit 443-681-3842 or 443-681-3793 To order collection kits and/or specimen collection supplies, contact:

Outfit Unit 443-681-3777 or Fax 443-681-3850

For Specific Test Requirements Refer to:
"Guide to Public Health Laboratory Services"
Available on line: dhmh.maryland.gov/laboratories

LABELING SPECIMENS/SAMPLES

Printed labels with all required patient information are recommended

Print patient name, date of birth Print date and time the specimen was collected

DO NOT cover expiration date of collection container

Write specimen source on collection container(s)

PACKAGING SPECIMENS FOR TRANSPORT

Never place specimens with different temperature requirements in the same biobag

Use one (1) biobag per temperature requirement

Review test request form to ensure all test(s) have been marked

Verify all specimens have been labeled

Place folded request form(s) in outer pouch of biobag

Multiple specimens from the same patient with the same temperature requirements must be packaged together in one (1) biobag

URINE SPECIMENS - Refrigerate PACKAGING AND SHIPPING

Double bag all urine specimens

Urine specimens require absorbent towel in biobag with specimen (express excess air before sealing)

Place bagged urine specimen in second biobag with all refrigerated specimens from the same patient

Place folded test request form(s) in outer pouch of second bag