



STATE OF MARYLAND

DHMH

Maryland Department of Health and Mental Hygiene

Larry Hogan, Governor - Boyd Rutherford, Lt. Governor - Dennis R. Schrader, Secretary

Laboratories Administration
Robert A. Myers, Ph.D., Director

DATE: January 9, 2017

TO: Medical Laboratory Directors, Local Health Officers, and Health Care Providers

FROM: Robert A. Myers, Ph.D.
Director, Laboratories Administration

RE: Updated Guidance and Instructions for Submission of Specimens for Suspected Zika Virus
Infection Testing at the Maryland DHMH Laboratory (January 2017)

On Monday, January 9, 2017, the Maryland DHMH Laboratory will begin using the CDC developed Trioplex multiplex real-time PCR (RT-PCR) assay to detect Zika virus RNA in clinical specimens that have been pre-authorized for testing. The Trioplex assay is approved for use by the federal Food and Drug Administration (FDA) under an Emergency Use Authorization (EUA) and has the added advantage of being able to simultaneously screen certain acute specimen types for both chikungunya and dengue virus RNA. The DHMH Laboratory has performed extensive verification and validation studies in preparation for this test transition and found that the Trioplex assay demonstrated similar performance characteristics when compared to the existing single-plex Zika, dengue and chikungunya RT-PCR assays that are currently in use in our laboratory.

The procedures for preapproval of Zika testing at the DHMH Laboratory remain unchanged (*see attached Guidelines and Instructions for Zika Testing*) and testing of at-risk patients will continue to be performed in accordance with the most current CDC and DHMH guidance and algorithm (<http://phpa.dhmh.maryland.gov/pages/zika.aspx>). However, the types and volumes of specimens used in the Trioplex assay will be somewhat different (*see attached Guidelines and Instructions for Zika Testing; Specimens Collection and Handling Table*). The Trioplex assay has been validated and approved for use with larger volumes of serum which improves the sensitivity of the assay. Additionally, whole blood (purple top tube: EDTA anticoagulant) and cerebrospinal fluid (CSF) can also be tested in the Trioplex assay. Some studies have indicated that testing of whole blood during the acute phase of Zika infections might improve the detection of viral RNA. **The use of plasma for Zika RT-PCR at the DHMH Laboratory is now discontinued.** Please be aware that a serum specimen must be submitted along with the other specimen types that can be tested in the Trioplex PCR assay. **Urine, CSF and whole blood specimens received without an accompanying serum cannot be tested in the Trioplex assay.**

The DHMH Laboratory will continue to test serum specimens using the CDC developed FDA EUA approved Zika IgM antibody capture enzyme-linked immunosorbent assay (MAC-ELISA) for the presence of IgM antibodies. Extensive cross-reactivity in flavivirus serological assays has been documented. Therefore, if specimens are reactive in the Zika MAC-ELISA an additional paired convalescent serum might be required to complete additional plaque reduction neutralization testing (PRNT) testing at the CDC Laboratory to possibly identify the most recently infecting flavivirus.

Contact the DHMH Molecular Diagnostic Laboratory at (443) 681-3924 or (443) 681-3905 during normal business hours (8:00AM to 4:30 PM weekdays) for questions regarding the transition to the Trioplex assay for Zika virus RNA testing and the Arbovirus Serology Laboratory at (443) 681-3937 or (443) 681-3932 for inquiries related to Zika virus IgM testing.

Howard Haft, MD, MMM, CPE, FACPE
David Blythe, MD, MPH

1770 Ashland Avenue • Baltimore, Maryland 21205
443-681-3800 • TTY for Disabled - Maryland Relay Service 1-800-735-2258
Toll Free 1-877-4MD-DHMH • Web Site: <http://maryland.gov/laboratories/>

PLEASE NOTE – REVISED CHART (April 2018)

Maryland Department of Health (MDH) Laboratories Administration Guidelines and Instructions for Zika Testing

Preapproval of Zika Test Results

An infectious disease consultation with MDH is still required to authorize specimens for Zika virus testing at the MDH Laboratory, prior to submitting specimens. Contact the MDH Infectious Disease Epidemiology and Outbreak Response Bureau at 410-767-6700 (or after hours at 410-795-7365) for consultation. Prior to contacting MDH, review of the current MDH Zika virus testing guidance at <http://phpa.health.maryland.gov/pages/zika.aspx> is highly recommended, as well as the CDC guidance found in the links below.

Specimen Collection and Handling

Testing Category	Specimen Type (see Notes 1 & 2)	Volume/Amount	Collect In:	Storage and Shipping Conditions (see Note 3 for storage > ± 7 days)
Symptomatic Adults and Children Asymptomatic Pregnant Women (If testing at commercial laboratory, check with that laboratory for specimen instructions.) Refer to: http://www.cdc.gov/zika/hc-providers/types-of-tests.html	Serum	3-5 ml (6-10 ml blood draw)	Red top, tiger top, or gold top serum separator tube	Refrigeration (2-8°C)
	Whole Blood	4-5 ml	Purple top EDTA tube	Refrigeration (2-8°C)
	Urine	5-10 ml	Leak proof, sterile urine cup; label as urine	Refrigeration (2-8°C)
	Cerebral Spinal Fluid (CSF)	1-2 ml	Leak proof, sterile tube or vial; label as CSF	Refrigeration (2-8°C)
Infants (within 2 days of birth) Refer to: CDC Guidelines on Collecting and Submitting Specimens at Time of Birth for Zika Virus Testing. http://www.cdc.gov/zika/hc-providers/test-specimens-at-time-of-birth.html	Serum	≥ 2 ml serum (≥ 4 ml blood draw)	Red top, tiger top, or gold top serum separator tube	Refrigeration (2-8°C)
	Urine	5-10 ml	Leak proof, sterile urine cup; label as urine	Refrigeration (2-8°C)
	Fixed Placenta, Fetal membrane, Umbilical Cord	1 inch square of: ♦ Umbilical cord ♦ Fetal membranes ♦ Placental disk edge ♦ Placental disk midsection ♦ Pathologic Lesions	Fix specimens in formalin; volume of formalin used should be as small as possible, but about 10x mass of tissue. (See Note 4 for further instructions.)	Room Temperature

Notes

¹A serum specimen must accompany all urine, CSF, or whole blood specimens or testing will not be performed.

²Plasma is not accepted for Zika testing at MDH.

³If specimens (except whole blood and fixed tissue) are to be held for longer than seven (7) days after collection until delivery to the testing lab, it is recommended to freeze to ≤ 20°C and ship frozen (on dry ice). Avoid repeating freezing and thawing cycles. Whole blood EDTA should not be frozen but refrigerated and tested within one week of collection. Fixed tissue should be held and shipped at room temperature.

⁴Place tissue collected according to the dimensions provided above in 10% buffered formalin for three days (72 hours). After fixation, if not paraffin-embedded, tissues SHOULD be transferred to 70% ethanol for long term storage and for shipping.

For detailed instructions on how to submit other specimen types (including amniotic fluid and semen) for Zika, dengue, chikungunya and other arboviral tests, contact the MD DHMH Laboratories at (443) 681-3924 or (443) 681-3937 during normal business hours from 8:00 a.m. - 4:30 p.m., Monday through Friday.

- a. **Name of Health Department Person Approving Testing:** Please record on the requisition the name of the DHMH or local health department person approving the testing.
- b. **Clinical Illness/Compatible clinical presentation:** e.g., rash, acute onset fever, conjunctivitis, arthralgia
- c. **Pertinent travel history:** Recent travel to a region where local transmission of Zika virus has been documented (an updated list is available at <http://www.cdc.gov/zika/geo/index.html>)
- d. **History of any previous flavivirus infection:** e.g., West Nile virus (WNV), dengue virus
- e. **Acute illness onset date:** contemporaneous with the travel exposures in areas of ongoing transmission (illness onset date ≤ 14 days after exposure)
- f. **Immunization history:** Yellow fever (YF), Japanese encephalitis (JE), or Tick-borne encephalitis (TE) vaccines

During normal business hours (8:30AM to 4:30 PM weekdays), for questions regarding the transition to the Triplex assay for Zika virus RNA testing, please contact the DHMH Molecular Diagnostic Laboratory at 443-681-3924 or 443- 681-3905 or the Arbovirus Serology Laboratory at 443-681-3937 or 443-681-3924 for inquiries related to Zika virus IgM testing.

Must complete the test request authorization information (This is where reports will be sent). Include the name of Healthcare provider who can legally order the test(s) in "Test Request Authorized by"

Request Arbovirus Travel-Associated Panel. Provide specimen source:

Indicate "S" for serum – (SST or aliquot) or whole clotted blood (red top)

Accompanying specimens:

Indicate "B" for whole unclotted blood with EDTA (Purple top) **UNSPUN**

Indicate "U" for urine. (Leak-proof sterile urine cup)

Indicate "CSF" for Cerebrospinal fluid (Leak-proof sterile tube or vial)

***Urine, Whole blood, and CSF MUST be submitted with an accompanying serum specimen.**

Complete patient's Travel history (location and dates), symptoms (or asymptomatic), pregnancy status (including weeks of gestation) vaccination history, & immune status

For questions on Zika Virus testing, please contact the lab: PCR: (443) 681-3923/3924 Serology: (443) 681-3932/3937

Laboratories Administration MD DHMH
1770 Ashland Ave. • Baltimore, MD 21205
443-681-3800 <http://dhmh.maryland.gov/laboratories/>
Robert A. Myers, Ph.D., Director

STATE LAB
Use Only



Patient's first & last names must be on the specimen container and exactly match the lab slip

SEROLOGICAL TESTING

Health Care Provider Address _____ City _____ County _____ State _____ Zip Code _____ Contact Name _____ Phone _____ Fax _____		Patient SSN (last 4 digits) _____ Last Name _____ First Name _____ Date of Birth (mm/dd/yyyy) _____ Address _____ City _____ State _____ Zip Code _____ County _____	
Test Requested Authorized by _____ Sex: ☐ Male ☐ Female ☐ Transgender Male ☐ Transgender Female Race: ☐ American Indian/Alaska Native ☐ Asian ☐ Black/African American ☐ Native Hawaiian/Other Pacific Islander ☐ White MRN/Case # _____ DOB # _____ Date Collected _____ Previous Test Date? ☐ no ☐ yes Name of Test _____ Name of Test _____ Onset Date _____ Exposure Date _____		*Specimen Source Code Arbovirus Panel (Serum or CSF) Mandatory: Onset Date, Collection Date, and Travel History Arbovirus Travel-Associated Panel (Chikungunya, Dengue) Based on information provided PCR and/or immunological assays will be performed. Required information, check all that apply: DIAGNOSIS: ☐ Dengue ☐ Chikungunya ☐ Zika Virus (Other) _____ SYMPTOMS: ☐ Headache ☐ Fever ☐ Rash ☐ Joint/muscle pain ☐ Fatigue ☐ Nausea/vomiting/diarrhea Unilateral mental status change weakness/Oral Other: _____ ALLNESS FATAL? ☐ yes ☐ no TRAVEL HISTORY (dates and places) _____ IMMUNIZATIONS: Yellow fever? ☐ yes ☐ no Polio? ☐ yes ☐ no Malaria? ☐ yes ☐ no Hepatitis A? ☐ yes ☐ no Hepatitis B? ☐ yes ☐ no Hepatitis C? ☐ yes ☐ no Hepatitis D? ☐ yes ☐ no Hepatitis E? ☐ yes ☐ no Rabies? ☐ yes ☐ no Tetanus? ☐ yes ☐ no Pertussis? ☐ yes ☐ no Measles/Mumps/Rubella? ☐ yes ☐ no Varicella? ☐ yes ☐ no Polio? ☐ yes ☐ no Hepatitis A? ☐ yes ☐ no Hepatitis B? ☐ yes ☐ no Hepatitis C? ☐ yes ☐ no Hepatitis D? ☐ yes ☐ no Hepatitis E? ☐ yes ☐ no	
*Specimen Source Code Arbovirus Panel (Serum or CSF) Mandatory: Onset Date, Collection Date, and Travel History Arbovirus Travel-Associated Panel (Chikungunya, Dengue) Based on information provided PCR and/or immunological assays will be performed. Required information, check all that apply: DIAGNOSIS: ☐ Dengue ☐ Chikungunya ☐ Zika Virus (Other) _____ SYMPTOMS: ☐ Headache ☐ Fever ☐ Rash ☐ Joint/muscle pain ☐ Fatigue ☐ Nausea/vomiting/diarrhea Unilateral mental status change weakness/Oral Other: _____ ALLNESS FATAL? ☐ yes ☐ no TRAVEL HISTORY (dates and places) _____ IMMUNIZATIONS: Yellow fever? ☐ yes ☐ no Polio? ☐ yes ☐ no Malaria? ☐ yes ☐ no Hepatitis A? ☐ yes ☐ no Hepatitis B? ☐ yes ☐ no Hepatitis C? ☐ yes ☐ no Hepatitis D? ☐ yes ☐ no Hepatitis E? ☐ yes ☐ no Rabies? ☐ yes ☐ no Tetanus? ☐ yes ☐ no Pertussis? ☐ yes ☐ no Measles/Mumps/Rubella? ☐ yes ☐ no Varicella? ☐ yes ☐ no Polio? ☐ yes ☐ no Hepatitis A? ☐ yes ☐ no Hepatitis B? ☐ yes ☐ no Hepatitis C? ☐ yes ☐ no Hepatitis D? ☐ yes ☐ no Hepatitis E? ☐ yes ☐ no		*Specimen Source Code Arbovirus Panel (Serum or CSF) Mandatory: Onset Date, Collection Date, and Travel History Arbovirus Travel-Associated Panel (Chikungunya, Dengue) Based on information provided PCR and/or immunological assays will be performed. Required information, check all that apply: DIAGNOSIS: ☐ Dengue ☐ Chikungunya ☐ Zika Virus (Other) _____ SYMPTOMS: ☐ Headache ☐ Fever ☐ Rash ☐ Joint/muscle pain ☐ Fatigue ☐ Nausea/vomiting/diarrhea Unilateral mental status change weakness/Oral Other: _____ ALLNESS FATAL? ☐ yes ☐ no TRAVEL HISTORY (dates and places) _____ IMMUNIZATIONS: Yellow fever? ☐ yes ☐ no Polio? ☐ yes ☐ no Malaria? ☐ yes ☐ no Hepatitis A? ☐ yes ☐ no Hepatitis B? ☐ yes ☐ no Hepatitis C? ☐ yes ☐ no Hepatitis D? ☐ yes ☐ no Hepatitis E? ☐ yes ☐ no Rabies? ☐ yes ☐ no Tetanus? ☐ yes ☐ no Pertussis? ☐ yes ☐ no Measles/Mumps/Rubella? ☐ yes ☐ no Varicella? ☐ yes ☐ no Polio? ☐ yes ☐ no Hepatitis A? ☐ yes ☐ no Hepatitis B? ☐ yes ☐ no Hepatitis C? ☐ yes ☐ no Hepatitis D? ☐ yes ☐ no Hepatitis E? ☐ yes ☐ no	

Collection Date and Onset of symptoms Date **MUST** be completed

If specimens other than whole blood, urine, serum, or CSF are being requested, please note type of specimen here, e.g.:

Fresh or Fixed Tissue
Amniotic Fluid

You must write "Zika Virus" to request testing

Include the name of the Local Health Department or DHMH Epidemiologist who approved testing

Zika Virus
Approved by: ####

Original

STATE LAB
Use Only

Laboratories Administration MD DHMH
1770 Ashland Ave. • Baltimore, MD 21205
443-681-3800 <http://dhmh.maryland.gov/laboratories/>
Robert A. Myers, Ph.D., Director



Department of Health
and Mental Hygiene

SEROLOGICAL TESTING

☐ EH ☐ FP ☐ MTY/PN ☐ NOD ☐ STD ☐ TB ☐ CD ☐ COR	
Health Care Provider	
Patient SS# (last 4 digits):	
Address	
Last Name ☐ SR ☐ JR ☐ Other	
City	County
First Name M.I.	
State	Zip Code
Date of Birth (mm/dd/yyyy) / /	
Contact Name:	
Address	
Phone#	Fax#
City	
County	
State	
Zip Code	
Test Request Authorized by:	
Sex: ☐ Male ☐ Female ☐ Transgender M to F ☐ Transgender F to M	
Ethnicity: Hispanic or Latino Origin? ☐ yes ☐ no	
Race: ☐ American Indian/Alaska Native ☐ Asian ☐ Black/African American ☐ Native Hawaiian/other Pacific Islander ☐ White	
MRN/Case #	DOC #
Outbreak #	Submitter Lab #
Date Collected:	Time Collected: ☐ am ☐ pm
*Vaccination History: _____	
Previous Test Done? ☐ no ☐ yes	
Name of Test	Date
☐ 1st ☐ 2nd ☐ 3rd	
State Lab Number: _____	
Name of Test	
Date	
☐ 1st ☐ 2nd ☐ 3rd	
State Lab Number: _____	
Onset Date: _____	
Exposure Date: _____	
☐ Clinical Illness/Symptoms	

↓ SPECIMEN SOURCE CODE Arbovirus Panels (Serum or CSF) Mandatory: Onset Date, Collection Date, and Travel History ☐ Arbovirus Endemic Panel (WNV, EEE, SLE, LAC) ☐ Arbovirus Travel-Associated Panel (Chikungunya, Dengue) Based on information provided PCR and/or immunological assays will be performed. Required information, check all that apply: DIAGNOSIS: ☐ Aseptic Meningitis ☐ Encephalitis ☐ Other _____ SYMPTOMS: ☐ headache ☐ fever ☐ stiff neck ☐ altered mental state ☐ muscle weakness ☐ rash ☐ Other _____ ILLNESS FATAL? ☐ yes ☐ no TRAVEL HISTORY (dates and places) IMMUNIZATIONS: Yellow fever? ☐ yes ☐ no Flavivirus? ☐ yes ☐ no IMMUNOCOMPROMISED? ☐ yes ☐ no	↓ SPECIMEN SOURCE CODE Herpes Simplex Virus (HSV) Types 1&2 Legionella Leptospira Lyme Disease *MMRV Immunity Screen: [Measles (Rubeola), Mumps, Rubella, Varicella (Chickenpox) IgG Ab only] Mononucleosis - Infectious *Mumps Immunity Screen Mycoplasma Rocky Mountain Spotted Fever (RMSF) *Rabies (RFFIT) (*List vaccination dates above) *Rubella Immunity Screen *Rubeola (Measles) Immunity Screen Schistosoma Strongyloides Syphilis - Previously treated? ☐ yes ☐ no Toxoplasma Tularemia Varicella Immunity Screen VDRL (CSF only) CDC/Other Test(s) Add'l Specimen Codes _____ Prior arrangements have been made with the following DHMH Labs Administration employee: _____ _____ _____ _____ Please Note Vaccination History above*	↓ SPECIMEN SOURCE CODE * LAVENDER TOP TUBE REQUIRED Hemoglobin Disorders Blood transfusion? (last 4 months) ☐ yes ☐ no Prenatal screen? ☐ yes ☐ no Father of baby screen? ☐ yes ☐ no Guardian's name if patient is a minor: _____ Name of mother of "at risk" baby: _____ RESTRICTED TEST Pre-Approved Submitters Only Submit a separate specimen for HIV Instructions go to: http://dhmh.maryland.gov/laboratories/ HIV Country of Origin _____ Rapid Test: ☐ Reactive ☐ Negative Date: _____ Specimen stored refrigerated (2°-8°c) after collection. ☐ yes ☐ no Specimen transported on cold packs ☐ yes ☐ no SPECIMEN SOURCE CODE: PLACE CODE IN BOX NEXT TO TEST B Blood (5 ml) CSF Cerebrospinal Fluid L Lavender Top Tube P Plasma S Serum (1 ml per test) UR Urine
--	---	--

Original