

Preanalytical Quality Assurance: Submitting Specimens to the State TB Laboratory

Richard Oatis
Supervisor, Mycobacteriology
State of Maryland Department of Health and Mental Hygiene
Labs Administration
(443) 681-3944
richard.oatis@maryland.gov

Preanalytical Quality Assurance

- Everything done before laboratory diagnostic testing to ensure the most accurate results possible.
- Our results can only be as good as the specimen and information we receive.
- Poor handling and storage of specimens, as well as improper specimen labeling and missing information on lab requisitions, can affect the accuracy of our results.


Filling out the lab requisition

- Submitter information, including ordering physician.
- Date of specimen collection.
- Source of specimen and tests requested.
- Patient information, as completely as possible.
- Make sure identifiers on lab slip match the identifiers on the specimen.

Infectious Agents Form

STATE LAB
Use Only

Laboratories Administration MDH
1770 Ashland Ave • Baltimore, MD 21205
443-681-3800 <http://health.maryland.gov/laborders/>
Robert A. Myers, Ph.D., Director



MARYLAND
Department of Health

INFECTIOUS AGENTS: CULTURE/DETECTION

TYPE OF PRINT REQUIRED INFORMATION OR PLACE LABELS ON BOTH COPIES	<input type="checkbox"/> BH <input type="checkbox"/> FP <input type="checkbox"/> MTY/PM <input type="checkbox"/> NOO <input type="checkbox"/> STD <input type="checkbox"/> TB <input type="checkbox"/> CD <input type="checkbox"/> COR			Patient SS # (last 4 digits):	
	Health Care Provider			Last Name <input type="checkbox"/> SR <input type="checkbox"/> JR <input type="checkbox"/> Other:	
	Address			First Name M.I.	
	City County			Date of Birth (mm/kl/yyyy) / /	
	State Zip Code			Address	
	Contact Name:			City County	
	Phone # Fax #			State Zip Code	
	Test Request Authorized by:				
	Sex: <input type="checkbox"/> Male <input type="checkbox"/> Female <input type="checkbox"/> Transgender M to F <input type="checkbox"/> Transgender F to M Ethnicity: Hispanic or Latino Origin? <input type="checkbox"/> Yes <input type="checkbox"/> No				
	Race: <input type="checkbox"/> American Indian/Alaska Native <input type="checkbox"/> Asian <input type="checkbox"/> Black/African American <input type="checkbox"/> Native Hawaiian/Other Pacific Islander <input type="checkbox"/> White				
MRN/Case #		DOC #	Outbreak #	Submitter Lab #	
Date Collected:		Time Collected:	<input type="checkbox"/> a.m. <input type="checkbox"/> p.m.	Onset Date: / /	
Reason for Test: <input type="checkbox"/> Screening <input type="checkbox"/> Diagnosis <input type="checkbox"/> Contact <input type="checkbox"/> Test of Cure <input type="checkbox"/> 2-3 Months Post Rx <input type="checkbox"/> Suspected Carrier <input type="checkbox"/> Isolate for ID <input type="checkbox"/> Release					
Therapy/Drug Treatment: <input type="checkbox"/> No <input type="checkbox"/> Yes Therapy/Drug Type: _____					
SPECIMEN SOURCE CODE		SPECIMEN SOURCE CODE		SPECIMEN SOURCE CODE	
BACTERIOLOGY		MYCOBACTERIOLOGY/AF/ITB		SPECIAL BACTERIOLOGY	
Bacterial Culture - Routine		AFB/ITB Culture and Smear		Legionella Culture	
Add'l Specimen Codes: _____		AFB/ITB Referred Isolate for ID		Leptospira	
Bordetella pertussis		M. tuberculosis referred isolate for genotyping		Mycoplasma (Outbreak Investigation Only)	
Group A Strep		Nuclear Acid Amplification Test for		RESTRICTED TESTS	
Group B Strep - Screen		M. tuberculosis Complex (GeneXpert)		Pre-approved submitters only	
C. difficile Toxin		PARASITOLOGY		Chlamydia trachomatis/GC NAAT	
Diphtheria		Blood Parasites: _____		**Norovirus (See comment on reverse)	
Foot/Oral Pathogens		Country visited outside US: _____		QuantIFERON	
IB (e.coli, C. perfringens, S. aureus)		Ova & Parasites		Incubation: Time began: _____ a.m./p.m.	
Gonorrhea Culture:		Immigrant? <input type="checkbox"/> Yes <input type="checkbox"/> No		Time ended: _____ a.m./p.m.	
Incubated? <input type="checkbox"/> Yes <input type="checkbox"/> No		Cryptosporidium		OTHER TESTS FOR	
Hours incubated: _____		Cyclospora/Isospora		INFECTIOUS AGENTS	
Add'l specimen Codes: _____		Microsporidium		Test Name: _____	
MRSA (nile out)		Pinworm		Prior arrangements have been made with the following MDH Labs Administration employee:	
VRE (nile out)		VIRUS/CHLAMYDIA			
ENTERIC INFECTIONS		Adenovirus?			
Campylobacter		Chlamydia trachomatis culture			
E. coli O157 typing/Shiga toxins		Cytomegalovirus (CMV)			
Enteric Culture - Routine		Enterovirus (Includes Echo & Coxsackie)			
(Salmonella, Shigella, E. coli O157, Campylobacter)		Herpes Simplex Virus (Types 1 & 2)			
Salmonella typing		Influenza (Types A & B) Rapid Flu Test:			
Shigella typing		Type: _____			
Vibrio		Result: <input type="checkbox"/> Negative <input type="checkbox"/> Positive			
Yersinia		Patient admitted to hospital? <input type="checkbox"/> Yes <input type="checkbox"/> No			
REFERENCE MICROBIOLOGY		Parainfluenza (Types 1, 2 & 3)?			
ABC's (BIDS) # _____		Respiratory Syncytial Virus (RSV)?			
Organism: _____		VARICELLA (VZV)			
Bacteria Referred Culture for ID		*MAY INCLUDE RESPIRATORY SCREENING PANEL			
Specify: _____		Comments: _____			
				SPECIMEN SOURCE CODES	
				PLACE CODE IN BOX NEXT TO TEST	
				B Blood SP Sputum	
				BW Bronchial Washing T Throat	
				CSF Cerebrospinal Fluid URE Urethra	
				CX Cerebrospinal Fluid UFV Urine (1 st void)	
				E Eye UCC Urine (Clean Catch)	
				F Feces V Vagina	
				N Nasopharynx/nassa W Wound	
				P Penis O Other: _____	
				R Rectum	

Submitter and patient information

TYPE OR PRINT REQUIRED INFORMATION OR PLACE LABELS ON BOTH COPIES	<input type="checkbox"/> EH <input type="checkbox"/> FP <input type="checkbox"/> MTY/PN <input type="checkbox"/> NOD <input type="checkbox"/> STD <input type="checkbox"/> TB <input type="checkbox"/> CD <input type="checkbox"/> COR		Patient SS # (last 4 digits):	
	Health Care Provider		Last Name <input type="checkbox"/> SR <input checked="" type="checkbox"/> JR <input type="checkbox"/> Other:	
	Address		First Name M.I.	
	City	County	Date of Birth (mm/dd/yyyy) / /	
	State	Zip Code	Address	
	Contact Name:		City County	
	Phone #	Fax #	State Zip Code	
	Test Request Authorized by:			
	Sex: <input type="checkbox"/> Male <input type="checkbox"/> Female <input type="checkbox"/> Transgender M to F <input type="checkbox"/> Transgender F to M		Ethnicity: Hispanic or Latino Origin? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	Race: <input type="checkbox"/> American Indian/Alaska Native <input type="checkbox"/> Asian <input type="checkbox"/> Black/African American <input type="checkbox"/> Native Hawaiian/Other Pacific Islander <input type="checkbox"/> White			
	MRN/Case #	DOC #	Outbreak #	Submitter Lab #
	Date Collected:	Time Collected:	<input type="checkbox"/> a.m. <input type="checkbox"/> p.m.	Onset Date: / /
	Reason for Test: <input type="checkbox"/> Screening <input type="checkbox"/> Diagnosis <input type="checkbox"/> Contact <input type="checkbox"/> Test of Cure <input type="checkbox"/> 2-3 Months Post Rx <input type="checkbox"/> Suspected Carrier <input type="checkbox"/> Isolate for ID <input type="checkbox"/> Release			
Therapy/Drug Treatment: <input type="checkbox"/> No <input type="checkbox"/> Yes		Therapy/Drug Type: _____ Therapy/Drug Date: ____/____/____		

Specimen type

SPECIMEN SOURCE CODES

PLACE CODE IN BOX NEXT TO TEST

B	Blood	SP	Sputum
BW	Bronchial Washing	T	Throat
CSF	Cerebrospinal Fluid	URE	Urethra
CX	Cervix/Endocervix	UFV	Urine (1 st Void)
E	Eye	UCC	Urine (Clean Catch)
F	Feces	V	Vagina
N	Nasopharynx/Nasal	W	Wound
P	Penis	O	Other: _____
R	Rectum		_____

Test requested

SPECIMEN SOURCE CODE	
	MYCOBACTERIOLOGY/AFB/TB
	AFB/TB Culture and Smear
	AFB/TB Referred Isolate for ID
	<i>M. tuberculosis</i> referred Isolate for genotyping
	Nuclear Acid Amplification Test for <i>M. tuberculosis</i> Complex (GeneXpert)

Sputum specimens

- First morning specimens are often best.
- 5 mL is the optimum volume.
- Watery specimens may not yield good results.

Other specimen collection considerations:

- Follow specimen collection instructions included with collection kit.
- Check that the patient has collected an appropriate specimen.
- Ensure that tube is closed tightly and properly and that biohazard “ziplock” bag is sealed properly.
- Refrigerate specimens as much as possible during storage/transport.
- Submit to lab ASAP – do not batch!

1 calendar day delivery time goal:

- “Benchmark is receipt within 1 day of specimen collection.”
- Expedite lab testing.
- Disallow for the growth of contaminating organisms.

Contaminated culture

TB Clinical

Microscopy Report

Fluorochrome -

AFB Not Found

Performed by: M. Plehn

Date: 1/11/16

Final Culture Report

Contaminated. Please submit another sample.

Performed by: A. Rivera

Date: 1/28/16

Ken Jost – Texas Department of State Health Services

transit time (days)	% of sputum specimens with AFB Culture result:				Total
	contam	neg	ntm	tb	
0	1%	59%	19%	22%	100%
1	2%	70%	8%	20%	100%
2	3%	70%	8%	19%	100%
3	2%	73%	8%	17%	100%
4	3%	74%	8%	15%	100%
5	3%	74%	7%	15%	100%
6	3%	75%	6%	15%	100%
7	5%	72%	8%	15%	100%
8	7%	70%	7%	16%	100%
9	5%	72%	8%	15%	100%
10	6%	74%	7%	14%	100%
>10	8%	71%	7%	14%	100%
NA	2%	74%	12%	11%	100%
Total	3%	73%	8%	17%	100%



Current turnaround times:

2017	2018 YTD	Description of turnaround times (TAT) for initial diagnostic specimens
		1. Promote rapid delivery of specimens. (TAT goal: Specimens should be received in the laboratory within 24 hours of specimen collection.) Report the percent of specimens received within 1, 2, and 3 calendar days.
37	29	% of specimens received within 1 calendar day.
58	56	% of specimens received within 2 calendar days.
79	74	% of specimens received within 3 calendar days.

How can we work together to improve?

- Check lab reports for date collected vs. date received.
- Identify changes to practices or workflow to expedite specimen submission.
- Notify the lab of any ways we may be able to assist.
- Be aware of specimens taking >7 days to arrive in the lab and take steps to prevent.

GeneXpert testing

- Only performed on sputum, bronch wash, trach aspirates, and tissue specimens.
- Liquid abscess drainage or aspirates and swabs not acceptable.
- Automatically performed on smear positive specimens from patients with no known history of TB.
- Testing of smear negative specimens must be requested verbally within 7 days of specimen receipt in the laboratory.

Drug Susceptibility Testing

- Performed automatically on all new isolates of TB.
- Repeated every 3 months if pan sensitive and every 2 months if there is any resistance.
- Second line drugs set up automatically if any resistance to first line drugs is detected.
- Results should be available about 2-3 weeks after culture ID.

Quantiferon

- Proper specimen handling and incubation necessary for accurate results.
- Incubator temperatures must be monitored closely and documented.
- Be aware of tube expiration dates!
- Tubes must be filled to indicated level.


QFT request

RESTRICTED TESTS	
Pre-approved submitters only	
	<i>Chlamydia trachomatis</i> /GC NAAT
	**Norovirus (See comment on reverse)
	QuantiFERON
	Incubation: Time began: _____ a.m./p.m.
	Time ended: _____ a.m./p.m.

QFT labeling



QFT labeling (cont.)



https://health.maryland.gov/laboratories/Pages/Tuberculosis-(TB)-Laboratory.aspx

HOME

DIVISIONS

CERTIFICATES OF LICENSURE

LABORATORY ADVISORY COMMITTEE

Quick Links

- › [Lab Hours](#)
- › [Holiday Schedule](#)
- › [Influenza Webpage](#)
- › [Arbovirus Travel-Associated Panel \(Zika, Dengue, Chikungunya\)](#)
- › [ARLN Webpage](#)
- › [Biosafety Outreach for Sentinel Laboratories](#)
- › [Requests to Access or Update Records](#)
- › [About the Labs](#)

Laboratory Guides

- › [Guide to Public Health Laboratory Services](#)
- › [Guide to Environmental Laboratory Services](#)


Laboratory Testing & Order Forms

- › [Infectious Agents Culture Detection Form and Instructions \(Fillable Form\)](#)
- › [Serological Testing Form and Instructions \(Fillable Form\)](#)

TB/ Mycobacteriology Lab:
443-681-3942

Laboratory Scientist Supervisor:
Rich Oatis
443-681-3944

Laboratory Scientist Lead:
Bryan Burall
443-681-3950



The main focus of our work in the Mycobacteriology Laboratory is to provide testing to aid in the diagnosis and control of TB in the State of Maryland. As a secondary goal, our lab also acts as a reference lab, assisting area hospital labs in the identification of a wide range of mycobacterial species.

Testing offered Monday through Friday includes:

§ *Culture and Smear*: We process a wide range of clinical specimens for AFB smear and culture. All smears are read by Fluorochrome staining with positives confirmed by Kinyoun. Culture is performed in both MGIT broth and Lowenstein-Jensen slants and held for seven (7) weeks incubation.

Microbiology Labs:

[Microbiology Home](#)

[Bioterrorism](#)

[Clinical Microbiology](#)

[Enteric, Diarrheal & Foodborne Diseases](#)

[Gonorrhea\(GC\)](#)

[Tuberculosis \(TB\)](#)

Mycobacteriology Menu:

[TB Lab Announcements](#)

[TB/ Mycobacteriology FAQs](#)

[Sputum Collection Instructions](#)

[Therapeutic Drug Monitoring \(TDM\) Instructions](#)

[Specimen Collection QA](#)



MARYLAND
Department of Health

Questions?

- TB Lab – (443) 681-3942
- Rich – (443) 681-3944